

Product datasheet for **RC213754L2V**

SBEM (MUCL1) (NM_058173) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SBEM (MUCL1) (NM_058173) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SBEM
Synonyms:	SBEM
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_058173
ORF Size:	270 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213754).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_058173.1
RefSeq Size:	396 bp
RefSeq ORF:	273 bp
Locus ID:	118430
UniProt ID:	Q96DR8
Cytogenetics:	12q13.2
Protein Families:	Secreted Protein
MW:	8.9 kDa



[View online »](#)

Gene Summary:

May play a role as marker for the diagnosis of metastatic breast cancer.[UniProtKB/Swiss-Prot Function]