

Product datasheet for **RC213758L3V**

Lipocalin like 1 protein (LCNL1) (NM_207510) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Lipocalin like 1 protein (LCNL1) (NM_207510) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Lipocalin like 1 protein
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_207510
ORF Size:	492 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213758).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_207510.2
RefSeq Size:	3706 bp
RefSeq ORF:	495 bp
Locus ID:	401562
UniProt ID:	Q6ZST4
Cytogenetics:	9q34.3
MW:	17.6 kDa



[View online »](#)