

Product datasheet for **RC213813A1V**

Human P2X5 (P2RX5) (NM_175080) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human P2X5 (P2RX5) (NM_175080) AAV Particle
Tag:	Myc-DDK
Symbol:	P2X5
Synonyms:	LRH-1; P2X5; P2X5R
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC213813 representing NM_175080 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGGGCAGGCGGGCTGCAAGGGGCTCTGCCTGTCGCTGTTGCGACTACAAGACCGAGAAGTATGTCATCG
CCAAGAACAAGAAGGTGGGCCTGCTGTACCGGCTGCTGCAGGCCTCCATCCTGGCGTACCTGGTCGTATG
GGTGTTCCTGATAAAGAAGGGTTACCAAGACGTCGACACCTCCCTGCAGAGTGCTGTATCACCAAAGTC
AAGGGCGTGGCCTTACCAACACCTCGGATCTTGGGCAGCGGATCTGGGATGTCGCCGACTACGTCATTC
CAGCCCAGAATGAAGGCATTCTGATGGCGCGTGCTCCAAGGACAGCGACTGCCACGCTGGGGAAGCGGT
TACAGCTGGAACCGGAGTGAAGACCGGCCGCTGCCTGCGGAGAGAGAATTGGCCAGGGGCACCTGTGAG
ATCTTTGCCTGGTGCCCGTTGGAGACAAGCTCCAGGCCGAGGAGCCATTCTGAAGGAGCCGAAGACT
TCACCATTTTCATAAAGAACCACATCCGTTTCCCAAATTCAACTTCTCCAACAATGTGATGGACGTCAA
GGACAGATCTTTCCTGAAATCATGCCACTTTGGCCCCAAGAACCCTACTGCCCATCTTCCGACTGGGC
TCCGTGATCCGCTGGGCCGGGAGCGACTTCAGGATATAGCCCTGGAGGGTGGCGTGATAGGAATTAATA
TTGAATGGAAGTGTGATCTTGATAAAGCTGCCTCTGAGTGCCACCCTCACTATTCTTTTAGCCGTCTGGA
CAATAAACTTTCAAAGTCTGTCTCCTCCGGGTACAATTGAGATTGCCAGATATTACCGAGACGCGACC
GGGGTGGAGTTCCGCACCCTGATGAAAGCCTACGGGATCCGCTTTGACGTGATGGTGAACGGCAAGGGTG
CTTTCTTCTGCGACCTGGTACTCATCTACCTCATCAAAAAGAGAGAGTTTACCGTGACAAGAAGTACGA
GGAAGTGAAGGGCCTAGAAGACAGTTCCAGGAGGCCGAGGACGAGGCATCGGGGCTGGGGCTATCTGAG
CAGCTCACATCTGGGCCAGGGCTGCTGGGGATGCCGGAGCAGCAGGAGCTGCAGGAGCCACCCGAGGCGA
AGCGTGGAAGCAGCAGTCAGAAGGGGAACGGATCTGTGTGCCACAGCTCCTGGAGCCCCACAGGAGCAC
G

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC213813 representing NM_175080 Red=Cloning site Green=Tags(s) MGQAGCKGLCLSLFDYKTEKYVIAKNKKVGLLYRLLQASILAYLVVWVFLIKKGYQDVTSLQSAVITKV KGVAFTNTSDLGQRIWDVADYVIPAQNEGIPDGACSKSDCHAGEAVTAGNGVKTGRCLRRNLARGTCE IFAWCPLETSSRPEEPFLKEAEDFTIFIKNHIRFPKFNFNNVMDVKDRSFLKSCHFGPKNHYCPIFRLG SVIRWAGSDFQDIALEGGVIGINIEWNCDLKAASECHPHYSFSRLDNKLSKSVSSGYNFRFARYYRDAA GVEFRTLKAYGIRFDVMVNGKAFFCDLVLIYLIKREFYRDKKYEEVRGLEDSSQEADEASGLGLSE QLTSGPGLLGMPQQELQEPPEAKRGSSSQKNGSVCPQLLEPHRST TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_175080
ORF Size:	1191 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_175080.1
RefSeq Size:	2150 bp
RefSeq ORF:	1194 bp
Locus ID:	5026
UniProt ID:	Q93086
Cytogenetics:	17p13.2
MW:	44.2 kDa