

Product datasheet for **RC213826L4V**

CRLS1 (NM_019095) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CRLS1 (NM_019095) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CRLS1
Synonyms:	C20orf155; CLS; CLS1; dj967N21.6; GCD10
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_019095
ORF Size:	903 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213826).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_019095.3
RefSeq Size:	1162 bp
RefSeq ORF:	906 bp
Locus ID:	54675
UniProt ID:	Q9UJA2
Cytogenetics:	20p12.3
Domains:	CDP-OH_P_transf
Protein Families:	Transmembrane



[View online »](#)

Protein Pathways: Glycerophospholipid metabolism, Metabolic pathways

MW: 32.4 kDa

Gene Summary: This gene encodes a member of the CDP-alcohol phosphatidyltransferase class-I family of proteins. The encoded enzyme catalyzes the synthesis of cardiolipin, a phospholipid component of mitochondrial membranes that is critical for mitochondrial function. [provided by RefSeq, Apr 2016]