

## Product datasheet for **RC213832L2V**

### **KBTBD5 (KLHL40) (NM\_152393) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	KBTBD5 (KLHL40) (NM_152393) Human Tagged ORF Clone Lentiviral Particle
Symbol:	KBTBD5
Synonyms:	KBTBD5; NEM8; SRYP; SYRP
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_152393
ORF Size:	1863 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213832).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_152393.2</a>
RefSeq Size:	2423 bp
RefSeq ORF:	1866 bp
Locus ID:	131377
UniProt ID:	<a href="#">Q2TBA0</a>
Cytogenetics:	3p22.1
Domains:	BTB, Kelch
MW:	69.1 kDa


[View online »](#)

**Gene Summary:**

This gene encodes a protein containing a BACK domain, a BTB/POZ domain, and 5 Kelch repeats, however, its exact function is not known. The gene and the multi-domain protein structure are conserved across different taxa, including primates, rodents, chicken and zebrafish. [provided by RefSeq, Dec 2012]