

Product datasheet for RC213837L3

PGAP1 (NM_024989) Human Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PGAP1 (NM_024989) Human Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	PGAP1
Synonyms:	Bst1; ISPD3024; MRT42; NEDDSBA; SPG67
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213837).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_024989
ORF Size:	2766 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_024989.3
RefSeq Size:	11113 bp
RefSeq ORF:	2769 bp
Locus ID:	80055
UniProt ID:	Q75T13
Cytogenetics:	2q33.1
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Glycosylphosphatidylinositol(GPI)-anchor biosynthesis, Metabolic pathways
MW:	105.4 kDa
Gene Summary:	The protein encoded by this gene functions early in the glycosylphosphatidylinositol (GPI) biosynthetic pathway, catalyzing the inositol deacylation of GPI. The encoded protein is required for the production of GPI that can attach to proteins, and this may be an important factor in the transport of GPI-anchored proteins from the endoplasmic reticulum to the Golgi. Defects in this gene are a cause an autosomal recessive form of cognitive impairment. [provided by RefSeq, Jul 2017]

Product images:



Circular map for RC213837L3