

## Product datasheet for **RC213837L4V**

### PGAP1 (NM\_024989) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	PGAP1 (NM_024989) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PGAP1
Synonyms:	Bst1; ISPD3024; MRT42; NEDDSBA; SPG67
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_024989
ORF Size:	2766 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213837).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_024989.3</a>
RefSeq Size:	11113 bp
RefSeq ORF:	2769 bp
Locus ID:	80055
UniProt ID:	<a href="#">Q75T13</a>
Cytogenetics:	2q33.1
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Glycosylphosphatidylinositol(GPI)-anchor biosynthesis, Metabolic pathways



[View online »](#)

**MW:** 105.4 kDa

**Gene Summary:** The protein encoded by this gene functions early in the glycosylphosphatidylinositol (GPI) biosynthetic pathway, catalyzing the inositol deacylation of GPI. The encoded protein is required for the production of GPI that can attach to proteins, and this may be an important factor in the transport of GPI-anchored proteins from the endoplasmic reticulum to the Golgi. Defects in this gene are a cause an autosomal recessive form of cognitive impairment. [provided by RefSeq, Jul 2017]