

## Product datasheet for **RC213846L2V**

### **KCNK4 (NM\_033310) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	KCNK4 (NM_033310) Human Tagged ORF Clone Lentiviral Particle
Symbol:	KCNK4
Synonyms:	FHEIG; K2p4.1; TRAAK; TRAAK1
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_033310
ORF Size:	1296 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213846).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_033310.1</a>
RefSeq Size:	1702 bp
RefSeq ORF:	1182 bp
Locus ID:	50801
UniProt ID:	<a href="#">Q9NYG8</a>
Cytogenetics:	11q13.1
Protein Families:	Druggable Genome, Ion Channels: Potassium, Transmembrane
MW:	42.5 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a member of the TWIK-related arachidonic acid-stimulated two pore potassium channel subfamily. The encoded protein homodimerizes and functions as an outwardly rectifying channel. This channel is regulated by polyunsaturated fatty acids, temperature and mechanical deformation of the lipid membrane. This protein is expressed primarily in neural tissues and may be involved in regulating the noxious input threshold in dorsal root ganglia neurons. Alternate splicing results in multiple transcript variants. Naturally occurring read-through transcripts also exist between this gene and the downstream testis expressed 40 (TEX40) gene, as represented in GenID: 106780802. [provided by RefSeq, Nov 2015]