

## Product datasheet for **RC213847L4V**

### PTP rho (PTPRT) (NM\_133170) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	PTP rho (PTPRT) (NM_133170) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PTP rho
Synonyms:	R-PTP-T; RPTP-rho; RPTPrho
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_133170
ORF Size:	2801 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213847).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_133170.2</a>
RefSeq Size:	12710 bp
RefSeq ORF:	4383 bp
Locus ID:	11122
UniProt ID:	<a href="#">O14522</a>
Cytogenetics:	20q12-q13.11
MW:	161.4 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene is a member of the protein tyrosine phosphatase (PTP) family. PTPs are known to be signaling molecules that regulate a variety of cellular processes including cell growth, differentiation, mitotic cycle, and oncogenic transformation. This PTP possesses an extracellular region, a single transmembrane region, and two tandem intracellular catalytic domains, and thus represents a receptor-type PTP. The extracellular region contains a meprin-A5 antigen-PTP (MAM) domain, Ig-like and fibronectin type III-like repeats. The protein domain structure and the expression pattern of the mouse counterpart of this PTP suggest its roles in both signal transduction and cellular adhesion in the central nervous system. Two alternatively spliced transcript variants of this gene, which encode distinct proteins, have been reported. [provided by RefSeq, Jul 2008]