

## Product datasheet for **RC213869L1V**

### USP54 (NM\_152586) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

|                           |                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|---------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:             | Lentiviral Particles                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Product Name:             | USP54 (NM_152586) Human Tagged ORF Clone Lentiviral Particle                                                                                                                                                                                                                                                                                                                                                                                       |
| Symbol:                   | USP54                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Synonyms:                 | bA137L10.3; bA137L10.4; C10orf29                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Mammalian Cell Selection: | None                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Vector:                   | pLenti-C-Myc-DDK (PS100064)                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Tag:                      | Myc-DDK                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| ACCN:                     | NM_152586                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| ORF Size:                 | 5052 bp                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(RC213869).                                                                                                                                                                                                                                                                                                                                                                                     |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.                                                                                                                                                                                                                                                                                                               |
| RefSeq:                   | <a href="#">NM_152586.3</a>                                                                                                                                                                                                                                                                                                                                                                                                                        |
| RefSeq Size:              | 6202 bp                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| RefSeq ORF:               | 5055 bp                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Locus ID:                 | 159195                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| UniProt ID:               | <a href="#">Q70EL1</a>                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Cytogenetics:             | 10q22.2                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Protein Families:         | Protease                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| MW:                       | 187.2 kDa                                                                                                                                                                                                                                                                                                                                                                                                                                          |



[View online »](#)

Gene Summary: Has no peptidase activity.[UniProtKB/Swiss-Prot Function]