

Product datasheet for **RC213907L3V**

ETV2 (NM_014209) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ETV2 (NM_014209) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ETV2
Synonyms:	ER71; ETSRP71
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_014209
ORF Size:	1026 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213907).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_014209.2
RefSeq Size:	1490 bp
RefSeq ORF:	1029 bp
Locus ID:	2116
UniProt ID:	O00321
Cytogenetics:	19q13.12
Protein Families:	ES Cell Differentiation/IPS
MW:	37.1 kDa



[View online »](#)

Gene Summary:

Binds to DNA sequences containing the consensus pentanucleotide 5'-CGGA[AT]-3'.
[UniProtKB/Swiss-Prot Function]