

## Product datasheet for **RC213949L4V**

### WNT2B (NM\_024494) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	WNT2B (NM_024494) Human Tagged ORF Clone Lentiviral Particle
Symbol:	WNT2B
Synonyms:	WNT13
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_024494
ORF Size:	1173 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213949).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_024494.1</a> , <a href="#">NP_078613.1</a>
RefSeq Size:	1970 bp
RefSeq ORF:	1176 bp
Locus ID:	7482
UniProt ID:	<a href="#">Q93097</a>
Cytogenetics:	1p13.2
Protein Families:	Secreted Protein



[View online »](#)

<b>Protein Pathways:</b>	Basal cell carcinoma, Hedgehog signaling pathway, Melanogenesis, Pathways in cancer, Wnt signaling pathway
<b>MW:</b>	43.6 kDa
<b>Gene Summary:</b>	This gene encodes a member of the wntless-type MMTV integration site (WNT) family of highly conserved, secreted signaling factors. WNT family members function in a variety of developmental processes including regulation of cell growth and differentiation and are characterized by a WNT-core domain. This gene may play a role in human development as well as carcinogenesis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2014]