

Product datasheet for **RC213957A1V**

Human ADPRM (NM_020233) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human ADPRM (NM_020233) AAV Particle
Tag:	Myc-DDK
Symbol:	ADPRM
Synonyms:	C17orf48; MDS006; NBLA03831
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC213957 representing NM_020233 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGATGATAAACCCCAATCCTGAAGCCCTAAGTGACAGTTCAGAGCGTCTTTCTCCTTTGGCGTCATCG
CAGATGTTCAATTTGCAGACTTAGAAGATGGCTTTAATTTCCAAGGAACCGCGCGGATCTACAGACA
TAGTCTTCTTCACTTACAGGGTGCCATTGAAGACTGGAATAATGAAAGCAGCATGCCCTGTTGTGTCCTT
CAGCTTGGAGATATCATCGATGGATATAATGCACAGTATAATGCATCCAAAAAGTCCCTAGAATTGTTA
TGGACATGTTCAAGAGGCTTAAAGTTCCAGTTCATACATACATGGGGAACCATGAATTCATAACTTCAG
TAGAGAGTATTTAACACACTCTAACTTAACACTAAGTTTCTAGAAGATCAGATTGTACATCATCTGAG
ACCATGCCTTCAGAAGATTATTATGCTTATCATTTTGTACCATTCCTAAATTCCGGTTCATTTTACTTG
ATGCATATGACTTGAGTGTCTTGGCGTGGATCAGTCTTCTCCAAAATACGAGCAGTGTATGAAGATATT
GAGGGAGCACAATCCAAATACGGAAGTGAATAGTCTCAAGGACTTTCTGAGCCCCAGTTTGTCCAGTTT
AATGGAGGATTCAGCCAAGAACAGCTAACTGGTTGAATGAAGTGCTAACATTCTGACACAAACCAAG
AAAAGGTGGTGATTGTGAGCCATCTTCCCATTTACCCGGACGCCTCTGACAATGTGTGCTGGCCTGGAA
CTACAGAGATGCCCTGGCAGTCATTTGGTCTCATGAGTGTGTGGTGTGTTTCTTGTGTCACACCCAT
GATGGTGGCTACTCTGAGGATCCTTTTGGTGTATACCAGTCAACCTAGAAGGAGTTATTGAAACAGCTC
CAGACAGCCAAGCCTTTGGCACAGTTCATGTCTATCCTGACAAAATGATGTTGAAAGGGAGAGGCAGAGT
TCCAGATAGAATTATGAATTACAAGAAAGAAAGAGCCTTCCATTGT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC213957 representing NM_020233 Red=Cloning site Green=Tags(s) MDDKPNPEALSDSSERLFSFGVIADVQFADLEDGFNFQGTRRRYYRHSLLHLQGAIEDWNNNESSMPCCVL QLGDIIDGYNAQYNASKKSLELVMDMFKRLKVPVHHTWGNHEFYNFSREYLTHSKLNTKFLEDQIVHHPE TMPSEDIYAYHFVPFPKFRFILLDAYDLSVLGVDQSSPKYEQCMKILREHNPNTLNSPQGLSEPQFVQF NGGFSQEQLNWLNEVLTFSDTNQEKVVIYSHLPYPDASDNVCLAWNYRDALAVIWSHECVVCFAGHTH DGGYSEDPFGVYHVNLEGVIIETAPDSQAFGTVHVYPDKMMLKGRGRVPDRIMNYKKERAFHC TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_020233
ORF Size:	1026 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_020233.4
RefSeq Size:	1518 bp
RefSeq ORF:	1029 bp
Locus ID:	56985
UniProt ID:	Q3LIE5
Cytogenetics:	17p13.1
MW:	39.3 kDa