

Product datasheet for **RC213963L4V**

EPHA8 (NM_001006943) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | EPHA8 (NM_001006943) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | EPHA8 |
| Synonyms: | E EK; EK3; HEK3 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_001006943 |
| ORF Size: | 1485 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC213963). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001006943.1 |
| RefSeq Size: | 1866 bp |
| RefSeq ORF: | 1488 bp |
| Locus ID: | 2046 |
| UniProt ID: | P29322 |
| Cytogenetics: | 1p36.12 |
| Protein Families: | Druggable Genome, Protein Kinase, Transmembrane |
| Protein Pathways: | Axon guidance |



[View online »](#)

MW: 53.9 kDa

Gene Summary: This gene encodes a member of the ephrin receptor subfamily of the protein-tyrosine kinase family. EPH and EPH-related receptors have been implicated in mediating developmental events, particularly in the nervous system. Receptors in the EPH subfamily typically have a single kinase domain and an extracellular region containing a Cys-rich domain and 2 fibronectin type III repeats. The ephrin receptors are divided into 2 groups based on the similarity of their extracellular domain sequences and their affinities for binding ephrin-A and ephrin-B ligands. The protein encoded by this gene functions as a receptor for ephrin A2, A3 and A5 and plays a role in short-range contact-mediated axonal guidance during development of the mammalian nervous system. [provided by RefSeq, Jul 2008]