

Product datasheet for **RC213969L3V**

FABP6 (NM_001040442) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	FABP6 (NM_001040442) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FABP6
Synonyms:	I-15P; I-BABP; I-BALB; I-BAP; ILBP; ILBP3; ILLBP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001040442
ORF Size:	531 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213969).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001040442.1 , NP_001035532.1
RefSeq Size:	672 bp
RefSeq ORF:	534 bp
Locus ID:	2172
UniProt ID:	P51161
Cytogenetics:	5q33.3
Protein Pathways:	PPAR signaling pathway
MW:	19.7 kDa



[View online »](#)

Gene Summary:

This gene encodes the ileal fatty acid binding protein. Fatty acid binding proteins are a family of small, highly conserved, cytoplasmic proteins that bind long-chain fatty acids and other hydrophobic ligands. FABP6 and FABP1 (the liver fatty acid binding protein) are also able to bind bile acids. It is thought that FABPs roles include fatty acid uptake, transport, and metabolism. Transcript variants generated by alternate transcription promoters and/or alternate splicing have been found for this gene. [provided by RefSeq, Jul 2008]