

Product datasheet for **RC213983A1V**

Human **OR2A14** (NM_001001659) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human OR2A14 (NM_001001659) AAV Particle
Tag:	Myc-DDK
Symbol:	OR2A14
Synonyms:	OR2A6; OR2A14P; OST182
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC213983 representing NM_001001659 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGAAGGCAACAAGACATGGATCACAGACATCACCTTGCCGCGATTCCAGGTTGGTCCAGCACTGGAGA
TTCTCCTCTGTGGACTTTTCTCTGCCTTCTATACACTCACCTGCTGGGAATGGGGTCATCTTTGGGAT
TATCTGCCTGGACTGTAAGCTTCACACCCCATGTACTTCTCCTCTCACACCTGGCCATTGTTGACATA
TCCTATGCTTCCAACATATGTCCCAAGATGCTGACGAATCTTATGAACCAGGAAAGCACCATCTCCTTTT
TTCCATGCATAATGCAGACATTCTGTATTTGGCTTTTGCTCACGTAGAGTGTCTGATTTTGGTGGTGAT
GTCCTATGATCGCTATGCGGACATCTGCCACCCCTTACGTTACAATAGCCTCATGAGCTGGAGAGTGTGC
ACTGTCCTGGCTGTGGCTTCTGGGTGTTGAGCTTCTCCTGGCTCTGGTCCCTTTAGTTCTCATCCTGA
GCCTGCCCTTCTGCGGGCCTCATGAAATCAACCACTTCTTCTGTGAAATCCTGTCTGCTCCTCAAGTTGGC
CTGTGCTGACACCTGGCTCAACCAGGTGGTCATCTTGCAGCCTGCGTGTTCATCCTGGTGGGGCCACTC
TGCTGGTGCTGGTCTCCTACTTGCGCATCTTGCCGCCCATCTTGAGGATCCAGTCTGGGAGGGCCGCA
GAAAGGCCTTCTCCACCTGCTCCTCCACCTTTGCGTGGTGGGACTCTTCTTTGGCAGCGCCATTGTCAC
GTACATGGCCCCAAGTCCCGCCATCTGAGGAGCAGCAGAAAGTTCTTCCCTGTTTTACAGCCTTTTC
AATCCAATGCTGAACCCCTGATATATAGCCTAAGGAATGCAGAGGTCAAGGGCGCCCTGAGGAGGGCAC
TGAGGAAGGAGAGGCTGACG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC213983 representing NM_001001659 Red=Cloning site Green=Tags(s) MEGNKTWITDITLPRFQVGPALIEILLCGLSAFYTLTLLGNGVIFGIICLDCKLHTPMYFFLSHLAIVDI SYASNYVPKMLTNLMNQESTISFFPCIMQTFLYLAFHVECLILVMSYDRYADICHPLRYNSLMSWRVC TVLAVASWVFSFLLALVPLVLILSLPFCGPHEINHFFCEILSVLKLACADTWLNQVVIFAACVFILVGPL CLVLVSYLRILAAILRIQSGEGRRKAFSTCSSHLCVVGLFFGSAIVTYMAPKSRHPPEEQQKVLSEFYSLF NPMLNPLIYSLRNAEVKGALRRALRKERLT TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001001659
ORF Size:	930 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001001659.1</u>
RefSeq Size:	933 bp
RefSeq ORF:	933 bp
Locus ID:	135941
UniProt ID:	<u>Q96R47</u>
Cytogenetics:	7q35
MW:	34.8 kDa