

## Product datasheet for **RC213996L2V**

### **GRAF (ARHGAP26) (NM\_015071) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	GRAF (ARHGAP26) (NM_015071) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GRAF
Synonyms:	GRAF; GRAF1; OPHN1L; OPHN1L1
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_015071
ORF Size:	2442 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213996).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_015071.2</a>
RefSeq Size:	6518 bp
RefSeq ORF:	2445 bp
Locus ID:	23092
UniProt ID:	<a href="#">Q9UNA1</a>
Cytogenetics:	5q31.3
Domains:	RhoGAP, SH3, PH
Protein Families:	Druggable Genome



[View online »](#)

MW: 92.1 kDa

**Gene Summary:** Interaction of a cell with the extracellular matrix triggers integrin cell surface receptors to begin signaling cascades that regulate the organization of the actin-cytoskeleton. One of the proteins involved in these cascades is focal adhesion kinase. The protein encoded by this gene is a GTPase activating protein that binds to focal adhesion kinase and mediates the activity of the GTP binding proteins RhoA and Cdc42. Defects in this gene are a cause of juvenile myelomonocytic leukemia (JMML). Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Mar 2017]