

Product datasheet for **RC213998L4V**

C2orf28 (ATRAID) (NM_080592) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	C2orf28 (ATRAID) (NM_080592) Human Tagged ORF Clone Lentiviral Particle
Symbol:	C2orf28
Synonyms:	APR--3; APR-3; APR3; C2orf28; HSPC013; p18; PRO240
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_080592
ORF Size:	852 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213998).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_080592.3 , NP_542159.3
RefSeq Size:	1278 bp
RefSeq ORF:	855 bp
Locus ID:	51374
UniProt ID:	Q6UW56
Cytogenetics:	2p23.3
Protein Families:	Druggable Genome, Transmembrane
MW:	30.3 kDa



[View online »](#)

Gene Summary:

This gene is thought to be involved in apoptosis, and may also be involved in hematopoietic development and differentiation. The use of alternative splice sites and promoters result in multiple transcript variants encoding different isoforms.[provided by RefSeq, Dec 2009]