

Product datasheet for **RC213999L1V**

KDM2B (NM_032590) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	KDM2B (NM_032590) Human Tagged ORF Clone Lentiviral Particle
Symbol:	KDM2B
Synonyms:	CXXC2; Fbl10; FBXL10; JHDM1B; PCCX2
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_032590
ORF Size:	4008 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC213999).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_032590.3
RefSeq Size:	5287 bp
RefSeq ORF:	4011 bp
Locus ID:	84678
UniProt ID:	Q8NHM5
Cytogenetics:	12q24.31
Protein Families:	Druggable Genome
MW:	152.4 kDa



[View online »](#)

Gene Summary:

This gene encodes a member of the F-box protein family which is characterized by an approximately 40 amino acid motif, the F-box. The F-box proteins constitute one of the four subunits of ubiquitin protein ligase complex called SCFs (SKP1-cullin-F-box), which function in phosphorylation-dependent ubiquitination. The F-box proteins are divided into 3 classes: Fbws containing WD-40 domains, Fbls containing leucine-rich repeats, and Fbxs containing either different protein-protein interaction modules or no recognizable motifs. The protein encoded by this gene belongs to the Fbls class. Multiple alternatively spliced transcript variants have been found for this gene, but the full-length nature of some variants has not been determined. [provided by RefSeq, Jul 2008]