

## Product datasheet for **RC214021L4V**

### **TMEM149 (IGFLR1) (NM\_024660) Human Tagged ORF Clone Lentiviral Particle**

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	TMEM149 (IGFLR1) (NM_024660) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TMEM149
Synonyms:	TMEM149
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_024660
ORF Size:	1065 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214021).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_024660.1</a>
RefSeq Size:	1217 bp
RefSeq ORF:	1068 bp
Locus ID:	79713
UniProt ID:	<a href="#">Q9H665</a>
Cytogenetics:	19q13.12
Protein Families:	Transmembrane
MW:	37.9 kDa


[View online »](#)

**Gene Summary:**

Probable cell membrane receptor for the IGF-like family proteins. Binds IGFL1 and IGFL3 with a higher affinity. May also bind IGFL2.[UniProtKB/Swiss-Prot Function]