

## Product datasheet for **RC214038L3V**

### EMILIN2 (NM\_032048) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	EMILIN2 (NM_032048) Human Tagged ORF Clone Lentiviral Particle
Symbol:	EMILIN2
Synonyms:	EMILIN-2; FOAP-10
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_032048
ORF Size:	3159 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214038).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_032048.2</a>
RefSeq Size:	4009 bp
RefSeq ORF:	3162 bp
Locus ID:	84034
UniProt ID:	<a href="#">Q9BXX0</a>
Cytogenetics:	18p11.32-p11.31
MW:	115.5 kDa



[View online »](#)

**Gene Summary:**

May be responsible for anchoring smooth muscle cells to elastic fibers, and may be involved not only in the formation of the elastic fiber, but also in the processes that regulate vessel assembly. Has cell adhesive capacity.[UniProtKB/Swiss-Prot Function]