

## Product datasheet for **RC214081L3V**

### DGKG (NM\_001346) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	DGKG (NM_001346) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DGKG
Synonyms:	DAGK3; DGK-GAMMA
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001346
ORF Size:	2373 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214081).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001346.1</a>
RefSeq Size:	3758 bp
RefSeq ORF:	2376 bp
Locus ID:	1608
UniProt ID:	<a href="#">P49619</a>
Cytogenetics:	3q27.2-q27.3
Domains:	DAGKa, DAGKc, EFh, DAG_PE-bind
Protein Families:	Druggable Genome



[View online »](#)

**Protein Pathways:** Glycerolipid metabolism, Glycerophospholipid metabolism, Metabolic pathways, Phosphatidylinositol signaling system

**MW:** 88.9 kDa

**Gene Summary:** This gene encodes an enzyme that is a member of the type I subfamily of diacylglycerol kinases, which are involved in lipid metabolism. These enzymes generate phosphatidic acid by catalyzing the phosphorylation of diacylglycerol, a fundamental lipid second messenger that activates numerous proteins, including protein kinase C isoforms, Ras guanyl nucleotide-releasing proteins and some transient receptor potential channels. Diacylglycerol kinase gamma has been implicated in cell cycle regulation and in the negative regulation of macrophage differentiation in leukemia cells. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]