

Product datasheet for **RC214084L4V**

NCR2 (NM_004828) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	NCR2 (NM_004828) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NCR2
Synonyms:	CD336; dj149M18.1; LY95; NK-p44; NKP44
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_004828
ORF Size:	828 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214084).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_004828.2 , NP_004819.1
RefSeq Size:	932 bp
RefSeq ORF:	831 bp
Locus ID:	9436
UniProt ID:	O95944
Cytogenetics:	6p21.1
Domains:	IG
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)

Protein Pathways: Natural killer cell mediated cytotoxicity

MW: 30.5 kDa

Gene Summary: Cytotoxicity-activating receptor that may contribute to the increased efficiency of activated natural killer (NK) cells to mediate tumor cell lysis.[UniProtKB/Swiss-Prot Function]