

Product datasheet for RC214085L4

MHF1 (CENPS) (NM_198544) Human Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MHF1 (CENPS) (NM_198544) Human Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	MHF1
Synonyms:	APITD1; APITD1-CORT; CENP-S; CENPS; FAAP16; MHF1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214085).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_198544
ORF Size:	492 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_198544.2
RefSeq Size:	1469 bp
RefSeq ORF:	495 bp
Locus ID:	100526739
UniProt ID:	Q8N2Z9
Cytogenetics:	1p36.22
MW:	18.9 kDa
Gene Summary:	This locus represents naturally occurring read-through transcription between the neighboring APITD1 (apoptosis-inducing, TAF9-like domain 1) and CORT (cortistatin) genes. Alternative splicing results in multiple transcript variants, two of which encode fusion proteins that share sequence identity with the products of each individual gene. [provided by RefSeq, Aug 2011]

Product images:



Circular map for RC214085L4