

Product datasheet for **RC214127L3V**

Muscarinic Acetylcholine Receptor M4 (CHRM4) (NM_000741) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Muscarinic Acetylcholine Receptor M4 (CHRM4) (NM_000741) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Muscarinic Acetylcholine Receptor M4
Synonyms:	HM4; M4R
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_000741
ORF Size:	1437 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214127).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_000741.2
RefSeq Size:	1468 bp
RefSeq ORF:	1440 bp
Locus ID:	1132
UniProt ID:	P08173
Cytogenetics:	11p11.2
Protein Families:	Druggable Genome



[View online »](#)

Protein Pathways: Neuroactive ligand-receptor interaction, Regulation of actin cytoskeleton

MW: 52.9 kDa

Gene Summary: The muscarinic cholinergic receptors belong to a larger family of G protein-coupled receptors. The functional diversity of these receptors is defined by the binding of acetylcholine and includes cellular responses such as adenylate cyclase inhibition, phosphoinositide degeneration, and potassium channel mediation. Muscarinic receptors influence many effects of acetylcholine in the central and peripheral nervous system. The clinical implications of this receptor are unknown; however, mouse studies link its function to adenylyl cyclase inhibition. [provided by RefSeq, Jul 2008]