

Product datasheet for **RC214137L4V**

MBNL2 (NM_144778) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	MBNL2 (NM_144778) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MBNL2
Synonyms:	MBLL; MBLL39; PRO2032
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_144778
ORF Size:	1101 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214137).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_144778.2
RefSeq Size:	4731 bp
RefSeq ORF:	1104 bp
Locus ID:	10150
UniProt ID:	Q5VZF2
Cytogenetics:	13q32.1
Domains:	zf-CCCH
Protein Families:	Transcription Factors


[View online »](#)

MW: 39.2 kDa

Gene Summary: This gene is a member of the muscleblind protein family which was initially described in *Drosophila melanogaster*. This gene encodes a C3H-type zinc finger protein that modulates alternative splicing of pre-mRNAs. Muscleblind proteins bind specifically to expanded dsCUG RNA but not to normal size CUG repeats and may thereby play a role in the pathophysiology of myotonic dystrophy. Several alternatively spliced transcript variants have been described but the full-length natures of only some have been determined. [provided by RefSeq, Mar 2012]