

Product datasheet for **RC214180A1V**

Human Cardiac Troponin T (TNNT2) (NM_000364) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human Cardiac Troponin T (TNNT2) (NM_000364) AAV Particle
Tag:	Myc-DDK
Symbol:	Cardiac Troponin T
Synonyms:	CMD1D; CMH2; CMPD2; cTnT; LVNC6; RCM3; TnTC
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC214180 representing NM_000364 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGCTCTGACATAGAAGAGGTGGTGAAGAGTACGAGGAGGAGGAGCAGGAAGAAGCAGCTGTTGAAGAAG
AGGAGGACTGGAGAGAGGACGAAGACGAGCAGGAGGAGGCAGCGGAAGAGGATGCTGAAGCAGAGGCTGA
GACCGAGGAGACCAGGGCAGAAGAAGATGAAGAAGAAGAGGAAGCAAAGGAGGCTGAAGATGGCCCAATG
GAGGAGTCCAAACCAAGCCAGGTCGTTTCATGCCCAACTTGGTGCCTCCCAAGATCCCGATGGAGAGA
GAGTGGACTTTGATGACATCCACCGGAAGCGCATGGAGAAGGACCTGAATGAGTTGCAGGCGCTGATTGA
GGCTCACTTTGAGAACAGGAAGAAAGAGGAGGAGGAGCTCGTTTCTCTCAAAGACAGGATCGAGAGACGT
CGGGCAGAGCGGGCCGAGCAGCAGCGCATCCGGAATGAGCGGGAGAAGGAGCGGCAGAACCGCCTGGCTG
AAGAGAGGGCTCGACGAGAGGAGGAGGAGAACAGGAGGAAGGCTGAGGATGAGGCCCCGAAGAAGAAGGC
TTTGTCCAACATGATGCATTTTGGGGTTACATCCAGAAGACAGAGCGGAAAAGTGGGAAGAGGCAGACT
GAGCGGAAAAGAAGAAGAAGATTCTGGCTGAGAGGAGGAAGGTGCTGGCCATTGACCACCTGAATGAAG
ATCAGCTGAGGGAGAAGGCCAAGGAGCTGTGGCAGAGCATCTATAACTTGGAGGCAGAGAAGTTCGACCT
GCAGGAGAAGTTCAAGCAGCAGAAATATGAGATCAATGTTCTCCGAAACAGGATCAACGATAACCCAGAAA
GTCTCCAAGACCCGCGGAAGGCTAAAGTCACCGGGCGCTGAAAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC214180 representing NM_000364 Red=Cloning site Green=Tags(s) MSDIEEVVEEYEEEEQEEAAVEEEEDWREDEDEQEEAAEEDAEAEATEETRAEEDEEEEEEAKEAEDGPM EESKPKPRSFMPNLVPPKIPDGERVDFDDIHRKRMEKDLNELQALIEAHFENRKKEEELVSLKDRIERR RAERAEQQRIRNEREKERQNLAEERARREEENRRKAEDARKKKALSNMMHFGGYIQKTERKSGKRQT EREKKKKILAERRKVLADHNLNEDQLREKAKELWQSIYNLEAEKFDLQEKFKQKYEINVLNRNRINDNQK VSKTRGKAKVTGRWK TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_000364
ORF Size:	885 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_000364.2
RefSeq Size:	1153 bp
RefSeq ORF:	888 bp
Locus ID:	7139
UniProt ID:	P45379
Cytogenetics:	1q32.1
MW:	35.4 kDa