

Product datasheet for **RC214186L4V**

DYX1C1 (DNAAF4) (NM_130810) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	DYX1C1 (DNAAF4) (NM_130810) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DNAAF4
Synonyms:	CILD25; DYX1; DYX1C1; DYXC1; EKN1; RD
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_130810
ORF Size:	1260 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214186).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_130810.1
RefSeq Size:	1993 bp
RefSeq ORF:	1263 bp
Locus ID:	161582
UniProt ID:	Q8WXU2
Cytogenetics:	15q21.3
Domains:	TPR
MW:	48.3 kDa



[View online »](#)

Gene Summary:

This gene encodes a tetratricopeptide repeat domain-containing protein. The encoded protein interacts with estrogen receptors and the heat shock proteins, Hsp70 and Hsp90. An homologous protein in rat has been shown to function in neuronal migration in the developing neocortex. A chromosomal translocation involving this gene is associated with a susceptibility to developmental dyslexia. Mutations in this gene are associated with deficits in reading and spelling. Alternative splicing results in multiple transcript variants. Read-through transcription also exists between this gene and the downstream cell cycle progression 1 (CCPG1) gene. [provided by RefSeq, Mar 2011]