

## Product datasheet for **RC214203L3V**

### ZC3H14 (NM\_207660) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ZC3H14 (NM_207660) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZC3H14
Synonyms:	MRT56; MSUT-2; NY-REN-37; SUT2; UKp68
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_207660
ORF Size:	1737 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214203).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_207660.2</a>
RefSeq Size:	3676 bp
RefSeq ORF:	1740 bp
Locus ID:	79882
UniProt ID:	<a href="#">Q6PJT7</a>
Cytogenetics:	14q31.3
MW:	65.2 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene is a poly(A)-binding protein that can affect gene expression and poly(A) tail length. The encoded protein may influence mRNA stability, nuclear export, and translation. [provided by RefSeq, May 2016]