

Product datasheet for **RC214207L1V**

NOX5 (NM_024505) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	NOX5 (NM_024505) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NOX5
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_024505
ORF Size:	2241 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214207).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_024505.2 , NP_078781.2
RefSeq Size:	2582 bp
RefSeq ORF:	2298 bp
Locus ID:	79400
UniProt ID:	Q96PH1
Cytogenetics:	15q23
Protein Families:	Ion Channels: Other, Transmembrane
MW:	84.5 kDa



[View online »](#)

Gene Summary:

This gene is predominantly expressed in the testis and lymphocyte-rich areas of spleen and lymph nodes. It encodes a calcium-dependent NADPH oxidase that generates superoxide, and functions as a calcium-dependent proton channel that may regulate redox-dependent processes in lymphocytes and spermatozoa. Alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Oct 2011]