

Product datasheet for **RC214211L4V**

EVX1 (NM_001989) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	EVX1 (NM_001989) Human Tagged ORF Clone Lentiviral Particle
Symbol:	EVX1
Synonyms:	EVX-1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001989
ORF Size:	1221 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214211).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001989.3
RefSeq Size:	1858 bp
RefSeq ORF:	1224 bp
Locus ID:	2128
UniProt ID:	P49640
Cytogenetics:	7p15.2
Protein Families:	ES Cell Differentiation/IPS
MW:	42.3 kDa



[View online »](#)

Gene Summary:

This gene encodes a member of the even-skipped homeobox family characterized by the presence of a homeodomain closely related to the *Drosophila* even-skipped (*eve*) segmentation gene of the pair-rule class. The encoded protein may play an important role as a transcriptional repressor during embryogenesis. [provided by RefSeq, Jul 2008]