

## Product datasheet for **RC214224L4V**

### **NARG1L (NAA16) (NM\_024561) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	NARG1L (NAA16) (NM_024561) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NARG1L
Synonyms:	NARG1L
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_024561
ORF Size:	2592 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214224).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_024561.3</a>
RefSeq Size:	4266 bp
RefSeq ORF:	2595 bp
Locus ID:	79612
UniProt ID:	<a href="#">Q6N069</a>
Cytogenetics:	13q14.11
Protein Families:	Druggable Genome, Transcription Factors
MW:	101.3 kDa



[View online »](#)

**Gene Summary:**

Auxillary subunit of the N-terminal acetyltransferase A (NatA) complex which displays alpha (N-terminal) acetyltransferase activity.[UniProtKB/Swiss-Prot Function]