

## Product datasheet for **RC214236L3V**

### LONRF3 (NM\_001031855) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	LONRF3 (NM_001031855) Human Tagged ORF Clone Lentiviral Particle
Symbol:	LONRF3
Synonyms:	RNF127
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001031855
ORF Size:	2277 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214236).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001031855.1</a>
RefSeq Size:	2625 bp
RefSeq ORF:	2280 bp
Locus ID:	79836
UniProt ID:	<a href="#">Q496Y0</a>
Cytogenetics:	Xq24
Protein Families:	Druggable Genome, Protease
MW:	84.3 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene contains a RING finger domain, a motif present in a variety of functionally distinct proteins and known to be involved in protein-protein and protein-DNA interactions. Multiple alternatively spliced transcript variants have been suggested, but their full length natures are not clear. [provided by RefSeq, Jul 2008]