

## Product datasheet for **RC214246**

### TSC22D3 (NM\_001015881) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** TSC22D3 (NM\_001015881) Human Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** TSC22D3  
**Synonyms:** DIP; DSIPI; GILZ; TSC-22R  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RC214246 representing NM\_001015881  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGAACACCGAAATGTATCAGACCCCCATGGAGGTGGCGGTCTACCAGCTGCACAATTTCTCCATCTCCT  
TCTTCTCTTCTGCTTGGAGGGGATGTGGTTCCGTTAAGCTGGACAACAGTGCCTCCGGAGCCAGCGT  
GGTGGCCATAGACAACAAGATCGAACAGGCCATGGATCTGGTGAAGAATCATCTGATGTATGCTGTGAGA  
GAGGAGGTGGAGATCCTGAAGGAGCAGATCCGAGAGCTGGTGGAGAAGAACTCCCAGCTAGAGCGTGAGA  
ACACCCTGTTGAAGACCCTGGCAAGCCCAGAGCAGCTGGAGAAGTTCCAGTCTGTCTGAGCCCTGAAGA  
GCCAGCTCCCGAATCCCCACAAGTGCCCGAGGCCCTGGTGGTTCTGCGGTG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC214246 representing NM\_001015881  
Red=Cloning site Green=Tags(s)

MNTEMYQTPMEVAVYQLHNFSISFFSSLLGGDVSVKLDNSASGASVVAIDNKIEQAMDLVKNHLMYAVR  
EEVEILKEQIRELVEKNSQLERENTLLKTLASPEQLEKFSQCLSPPEEPAPESPQVPEAPGGSVAV

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM\_001015881

ORF Size: 405 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001015881.1</a></u> , <u><a href="#">NP_001015881.1</a></u>
<b>RefSeq Size:</b>	1781 bp
<b>RefSeq ORF:</b>	234 bp
<b>Locus ID:</b>	1831
<b>UniProt ID:</b>	<u><a href="#">Q99576</a></u>
<b>Cytogenetics:</b>	Xq22.3
<b>Protein Families:</b>	Transcription Factors
<b>MW:</b>	14.8 kDa
<b>Gene Summary:</b>	This gene encodes the anti-inflammatory protein glucocorticoid (GC)-induced leucine zipper. Expression of this gene stimulated by glucocorticoids and interleukin 10 and it appears to play a key role in the anti-inflammatory and immunosuppressive effects of this steroid. This protein has also been shown to inhibit pro-inflammatory molecules including nuclear factor $\kappa$ B. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Jan 2016]