

## Product datasheet for **RC214265L1V**

### AP5B1 (NM\_138368) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	AP5B1 (NM_138368) Human Tagged ORF Clone Lentiviral Particle
Symbol:	AP5B1
Synonyms:	AP-5; PP1030
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_138368
ORF Size:	2634 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214265).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_138368.4</a>
RefSeq Size:	6595 bp
RefSeq ORF:	2637 bp
Locus ID:	91056
UniProt ID:	<a href="#">Q2VPB7</a>
Cytogenetics:	11q13.1
MW:	94.4 kDa



[View online »](#)

**Gene Summary:**

As part of AP-5, a probable fifth adaptor protein complex it may be involved in endosomal transport.[UniProtKB/Swiss-Prot Function]