

Product datasheet for **RC214266**

PI 3 Kinase p85 alpha (PIK3R1) (NM_181504) Human Tagged ORF Clone

Product data:

| | |
|---------------------------|---|
| Product Type: | Expression Plasmids |
| Product Name: | PI 3 Kinase p85 alpha (PIK3R1) (NM_181504) Human Tagged ORF Clone |
| Tag: | Myc-DDK |
| Symbol: | PI 3 Kinase p85 alpha |
| Synonyms: | AGM7; GRB1; IMD36; p85; p85-ALPHA |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |



[View online »](#)

ORF Nucleotide Sequence:

>RC214266 representing NM_181504
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGTACAATACTGTTTGGAAATATGGAAGACCTGGATTTAGAATATGCCAAGACAGATATAAATTGTGGCA
 CAGACTTGATGTTTTATATAGAAATGGACCCACCAGCACTGCCTCTAAACCACCAAAACCTACTACTGT
 AGCCAACAACGGTATGAATAACAATATGTCCTTACAAGATGCTGAATGGTACTGGGGAGATATCTCGAGG
 GAAGAAGTGAATGAAAACTTCGAGATACAGCAGACGGGACCTTTTTGGTACGAGATGCGTCTACTAAAA
 TGCATGGTGATTACTCTTACACTAAGGAAAGGGGAAATAACAAATTAATCAAATATTTTCATCGAGA
 TGGGAAATATGGCTTCTGACCCATTAACTTCAGTTCTGTGGTTGAATTAATAAACCCTACCGGAAT
 GAATCTCTAGCTCAGTATAATCCCAATGGATGTGAAATTAATCCAGTATCCAAATACCAACAGG
 ATCAAGTTGTCAAAGAAGATAATTGAAGCTGTAGGAAAAAATACATGAATATAACTCAGTTTCA
 AGAAAAAGTCGAGAATATGATAGATTATGAAGAATATACCCGCACATCCAGGAAATCCAAATGAAA
 AGGACAGCTATTGAAGCATTTAATGAAACCATAAAAAATTTGAAGAAGCAGTGCCAGACCAAGAGCGGT
 ACAGCAAAGAATACATAGAAAAGTTAAACGTGAAGGCAATGAGAAAGAATACAAGGATTATGCATAA
 TTATGATAAGTTGAAGTCTCGAATCAGTGAAATTAATGACAGTAGAAGAAGATTGGAAGAAGACTTGAAG
 AAGCAGGCAGCTGAGTATCGAGAAATGACAAACGTATGAACAGCATTAAACCAGACCTTATCCAGCTGA
 GAAAGACGAGAGACCAATACTTGATGTGGTTGACTCAAAAAGGTGTTCCGCAAAAGAAGTTGAACGAGTG
 GTTGGGCAATGAAACACTGAAGACCAATTAATCACTGGTGAAGATGATGAAGATTTGCCCATCATGAT
 GAGAAGACATGGAATGTTGGAAGCAGCAACCGAAACAAAGCTGAAAACCTGTTGCGAGGGGAAAGCGAGATG
 GCACTTTTCTTCCGGGAGAGCAGTAAACAGGGCTGCTATGCCTGCTGTAGTGGTGGACGGCGAAGT
 AAAGCATTGTGTCATAAAACAAAACAGCAACTGGCTATGGCTTTGCCGAGCCCTATAAATTTGTACAGCTCT
 CTGAAAGAAGTGGTGTACATTACCAACACACCTCCCTTGTGCAGCACAACGACTCCCTCAATGTACAC
 TAGCTACCCAGTATATGCACAGCAGAGGCGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC214266 representing NM_181504
 Red=Cloning site Green=Tags(s)

MYNTVWNMEDLDLEYAKTDINCGTDLMFYIEMDPPALPPKPPKPTTVANNGMNNMSLQDAEWYWGDISR
 EEVNEKLRDADGTFVLRDASTKMHGDYTLTRKGGNNKLIKIFHRDGKYGFSDPLTFSSVVELINHYRN
 ESLAQYNPKLDVKLLYPVSKYQQDQVVKEDNIEAVGKKLHEYNTQFQEKSRDYRLYEEYTRTSQEIQMK
 RTAIEAFNETIKIFEEQCQTQERYKEYIEKFKREGNEKEIQRIMHNYDKLSRISEIIDSRRRLEEDLK
 KQAAEYREIDKRMNSIKPDLIQLRKRTRDQYLMWLTQKGVRRQKLNWLNENTEDQYSLVEDDEDLPHHD
 EKTWNVGSNNRKAENLLRGKRDGTFVRESSKQGCYACSVVVDGEVKHCVINKTATGYGFAEPYNLYSS
 LKELVLHYQHTSLVQHNDLNVTLAYPVYAQRR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms:

https://cdn.origene.com/chromatograms/mk6122_a05.zip

Restriction Sites:

Sgfl-Mlul

Cloning Scheme:


ACCN: NM_181504

ORF Size: 1362 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_181504.4](#)

RefSeq Size: 5663 bp

RefSeq ORF: 1365 bp

Locus ID: 5295

UniProt ID: [P27986](#)

Cytogenetics: 5q13.1


Protein Families: Druggable Genome

Protein Pathways: Acute myeloid leukemia, Apoptosis, B cell receptor signaling pathway, Chemokine signaling pathway, Chronic myeloid leukemia, Colorectal cancer, Endometrial cancer, ErbB signaling pathway, Fc epsilon RI signaling pathway, Fc gamma R-mediated phagocytosis, Focal adhesion, Glioma, Insulin signaling pathway, Jak-STAT signaling pathway, Leukocyte transendothelial migration, Melanoma, mTOR signaling pathway, Natural killer cell mediated cytotoxicity, Neurotrophin signaling pathway, Non-small cell lung cancer, Pancreatic cancer, Pathways in cancer, Phosphatidylinositol signaling system, Progesterone-mediated oocyte maturation, Prostate cancer, Regulation of actin cytoskeleton, Renal cell carcinoma, Small cell lung cancer, T cell receptor signaling pathway, Toll-like receptor signaling pathway, Type II diabetes mellitus, VEGF signaling pathway


MW: 53.3 kDa

Gene Summary: Phosphatidylinositol 3-kinase phosphorylates the inositol ring of phosphatidylinositol at the 3-prime position. The enzyme comprises a 110 kD catalytic subunit and a regulatory subunit of either 85, 55, or 50 kD. This gene encodes the 85 kD regulatory subunit. Phosphatidylinositol 3-kinase plays an important role in the metabolic actions of insulin, and a mutation in this gene has been associated with insulin resistance. Alternative splicing of this gene results in four transcript variants encoding different isoforms. [provided by RefSeq, Jun 2011]


Product images:



Circular map for RC214266



Western blot validation of overexpression lysate (Cat# [LY403621]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC214266 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified PIK3R1 protein (Cat# [TP314266]). The protein was produced from HEK293T cells transfected with PIK3R1 cDNA clone (Cat# RC214266) using MegaTran 2.0 (Cat# [TT210002]).