

Product datasheet for **RC214292L1V**

CECR1 (ADA2) (NM_177405) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | CECR1 (ADA2) (NM_177405) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | ADA2 |
| Synonyms: | ADGF; CECR1; IDGFL; PAN; SNEDS; VAIHS |
| Mammalian Cell Selection: | None |
| Vector: | pLenti-C-Myc-DDK (PS100064) |
| Tag: | Myc-DDK |
| ACCN: | NM_177405 |
| ORF Size: | 810 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC214292). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_177405.1 |
| RefSeq Size: | 3071 bp |
| RefSeq ORF: | 813 bp |
| Locus ID: | 51816 |
| UniProt ID: | Q9NZK5 |
| Cytogenetics: | 22q11.1 |
| Protein Families: | Druggable Genome, Secreted Protein, Transmembrane |
| MW: | 30.5 kDa |



[View online »](#)

Gene Summary:

This gene encodes a member of a subfamily of the adenosine deaminase protein family. The encoded protein is one of two adenosine deaminases found in humans, which regulate levels of the signaling molecule, adenosine. The encoded protein is secreted from monocytes undergoing differentiation and may regulate cell proliferation and differentiation. This gene may be responsible for some of the phenotypic features associated with cat eye syndrome. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2013]