

Product datasheet for **RC214294L3V**

Exonuclease 1 (EXO1) (NM_006027) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Exonuclease 1 (EXO1) (NM_006027) Human Tagged ORF Clone Lentiviral Particle
Symbol:	EXO1
Synonyms:	HEX1; hExol
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_006027
ORF Size:	2538 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214294).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_006027.3 , NP_006018.3
RefSeq Size:	3210 bp
RefSeq ORF:	2541 bp
Locus ID:	9156
UniProt ID:	Q9UQ84
Cytogenetics:	1q43
Protein Families:	Druggable Genome, Stem cell - Pluripotency
Protein Pathways:	Mismatch repair




[View online »](#)

MW: 94 kDa

Gene Summary: This gene encodes a protein with 5' to 3' exonuclease activity as well as an RNase H activity. It is similar to the *Saccharomyces cerevisiae* protein Exo1 which interacts with Msh2 and which is involved in mismatch repair and recombination. Alternative splicing of this gene results in three transcript variants encoding two different isoforms. [provided by RefSeq, Jul 2008]

Product images:



[RC214294L3] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC214294L3V particle to overexpress human EXO1-Myc-DDK fusion protein.