

## Product datasheet for **RC214352L4V**

### CCNB1IP1 (NM\_021178) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	CCNB1IP1 (NM_021178) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CCNB1IP1
Synonyms:	C14orf18; HEI10
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_021178
ORF Size:	831 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214352).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_021178.2</a>
RefSeq Size:	1499 bp
RefSeq ORF:	834 bp
Locus ID:	57820
UniProt ID:	<a href="#">Q9NPC3</a>
Cytogenetics:	14q11.2
MW:	31.4 kDa



[View online »](#)

**Gene Summary:**

HEI10 is a member of the E3 ubiquitin ligase family and functions in progression of the cell cycle through G(2)/M.[supplied by OMIM, Apr 2004]