

## Product datasheet for **RC214357L1V**

### **GPR115 (ADGRF4) (NM\_153838) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	GPR115 (ADGRF4) (NM_153838) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GPR115
Synonyms:	GPR115; PGR18
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_153838
ORF Size:	2085 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214357).
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at <a href="mailto:custsupport@origene.com">custsupport@origene.com</a> or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a></p>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_153838.2</a> , <a href="#">NP_722580.2</a>
RefSeq Size:	2574 bp
RefSeq ORF:	2088 bp



[View online »](#)

<b>Locus ID:</b>	221393
<b>UniProt ID:</b>	<a href="#">Q8IZF3</a>
<b>Cytogenetics:</b>	6p12.3
<b>Protein Families:</b>	Druggable Genome, Transmembrane
<b>MW:</b>	77.5 kDa
<b>Gene Summary:</b>	Sequence analysis of this gene suggests that it encodes a member of the superfamily of G protein-coupled receptors. G protein-coupled receptors typically contain seven hydrophobic transmembrane domains, interact with guanine nucleotide binding regulatory proteins, and detect molecules outside the cell and act to transduce these signals into intracellular responses. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2016]