

Product datasheet for **RC214360**

PCDHGA6 (NM_018919) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PCDHGA6 (NM_018919) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	PCDHGA6
Synonyms:	PCDH-GAMMA-A6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC214360 representing NM_018919
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCGCCTCCGAGAGGCATCCGAGCGCAGCGAGCAGGTCCTGCTCCTCACGCTCCTGGGGACGCTGT
 GGGGGGCCGCGCAGCGCAGATCCGCTACTCTATTCCCGAGGAGCTGGAGAAAGGCTCCTTCGTAGGCAA
 CATCGTCAAGGATCTGGGACTGGAGCCCCAGGAGTTGGCGGAGCACGGAGTCCGCATCGTCTCCAGAGGT
 AGGATGCAGCTTTTCTCTCTGAATCCGCGAAACGGCAGCTTGGTCACCGCGGGTAGGATAGACCGCAGG
 AGCTCTGTGCTCAGAGCCCGGGTGTCTGGTGAGTTTAAACATCCTTGTGAGGATAAACTGAATCTTTA
 TCCCGTGGAAAGTGGAAATAGTGGACATTAATGACAATACACCCGATTCTTAAAGGAAGAATTGGAAGTG
 AAAATTCTCGAAACGCAGCTCCATCCTCTCGTTTTCCACTAATGGAGGTCTATGACCCTGATGTGGGAA
 TGAATCCCTTCAGGGATTTAAGCTCAGTGGTAATAGTCACTTCTCAGTGGACGTGCAAAGCGAAGCCCA
 TGGGCCCAAGTACCCGGAGCTGGTGTGGAGGGCACACTGGACCGGGAAGGAGAAGCCGTTTACCGCTG
 GTCTTACTGCCATGGATGGCGGCGACCTGTCCGCTCAAGCGTCGCCCAAATTCGTGTAACAGTTCTAG
 ATGTGAATGACAACACTCCAATGTTTACTCAGCCTGTCTACCGTGTAAAGTGTTCGAAAACCTGCCAGT
 AGGCACACCAGTGTGGCAGTACTGCCACCGACCAGGATGAAGGAGTCCACGGGGAAGTAACTTATTCC
 TTTGTGAAGATTACAGAAAAGATCTCACAATTTTCTGTTTGAATGTTTGTGAGGAAATTTCAACTT
 CTGCAATCTAGACTATGAGGACTCGAGTTTTATGAGCTGGGTGTTGAAGCCCGGGATGGGCCAGGTCT
 TCGAGACAGAGCGAAAGTCTTAATAACTATCTGGATGTCAATGATAATGTACCAGAAGTGGTTGTTACA
 TCTGGAAGCAGAACAAATGTCTGAAAGTGACCTCCAGGAACAGTAATCGCCCTTTTCAAGTGTTCGATC
 GAGACTCTGGCCTGAATGGCCTGGTAACCTGTTCCATCCCGAGAAGTCTCCCATTTGAATTGGAAAAATC
 AGTTGGCAATTATTATCGATTAGTGACAAATGCAGCTCTAGACCGGGAAGAGGTATTCTTGTAACATC
 ACTGTGACAGCCACGGACAAAGGAACACCACTGTCTACAGAAACAATCATCTCTCTAAATGTGGCAG
 ACACCAACGACAACCCGCCACCTTCCCCATTCTACTCAGTCTATGTCTTGAACAACCCAG
 GGGTGCCTCCATCTTCTGTGAATGACTGGACCCTGACGTGGACCAGAACGCCAAGTCTCTACTCA
 CTGGCAGAAGACACCTCCAGGGGGCGCCCTGTCTCTACGTGTCCATCAACTCCGACACTGGGATTC
 TGTACGCCCTCGCTCCTTCGACTATGAGCAGTTGAGAGACCTACAGCTGTGGGTGACAGCCAGCGACAG
 CGGGGACCCGCTCTTAGCAGCAACGTGTCACTGAGCCTGTTGTGCTGGACCAGAATGACAATGCGCCC
 GAGATCCTGTACCCGCCCTCCACAGACGGTCCACTGGCGTGGAGCTGGCGCCCGCTCCGAGAGC
 CCGGCTACCTGGTGACCAAGGTGGTGGCGGTGGACAGAGACTCGGGCCAGAACGCCTGGCTGTCTACCG
 CCTGCTCAAGGCCAGCGAGCCAGGACTTTTCTCAGTGGGCCTGCACACGGGCGAGGTGCGCACGGCGCGC
 GCCCTGTGGACAGAGACGCGCTCAAGCAGAGCCTAGTGGTGGCCGTCCAGGACCACGGCCAGCCCCCTC
 TCTCCGCCACTGTACGCTCACCGTGGCCGTGGCCGACAGGATCCCCGACATCCTGGCCGACCTGGGCGAG
 CCTCGAGCCCTCCGCCAAACCAACGATTCCGACCTCACTCTGTACCTGGTGGTGGCGGTGGCCGCGGT
 TCCTGCGTCTTCTGGCCTTCGTATCGTGTCTGGCGCTCAGACTGCAGCGCTGGCACAAGTCACGCC
 TGCTGCAGGCTTCGGGAGGTGGCTTAGCGAGCATGCCCGGCTCGCACTTTGTGGCGGTGGAAGGGTTTCG
 GGCTTTCTGCAGACCTATCCACGAGGTCTCACTCACTGCAGACTCGCGTAAGAGTCATCTGATTTTC
 CCCCAGCCCAACTATGCCGACACGCTTATCAACCAGGAGAGCTATGAGAAAAGCGAGCCTCTTCTGATAA
 CTCAGGATTTACTTGAACGAAAGGAGAACCCAGGCAACTTCAGCAAGCCCCGCCAACACGGACTGGCG
 TTTCTCTCAGGCCAGAGACCCGGCACCAGCGGCTCCAAAATGGCGATGACACCGGCACCTGGCCCAAC
 AACCAGTTTGACACAGAGATGCTGCAAGCCATGATCTTGGCGTCCGCCAGTGAAGCTGCTGATGGGAGCT
 CCACCCTGGGAGGGGTGCCGGCACCATGGGATTGAGCGCCCGCTACGGACCCAGTTACCCTGCAGCA
 CGTGCCCGACTACCGCCAGAATGTCTACATCCAGGCAGCAATGCCCACTGACCAACGCAGCTGGCAAG
 CGGGATGGCAAGGCCCCAGCAGGTGGCAATGGCAACAAGAAGAAGTCGGGCAAGAAGGAGAAGAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC214360 representing NM_018919
 Red=Cloning site Green=Tags(s)


MAPPQRHPQRSEQVLLLTLGLTWGAAAQIRYSIPEELEKGSFVGNIVKDLGLEPQELAEHGVRIVSRG
 RMQLFSLNPRNGSLVTAGRIDREELCAQSPRCLVSFNILVEDKLNLYPVEVEIVDINDNTPRFLKEELEV
 KILENAAPSSRFPLMEVYDPPDVGMSLQGFKLSGNSHFSVDVQSEAHGPKYPELVLEGLDREGEAVYRL
 VLTAMDGGDPVRSSVAQILVTVLVDVNDNTPMFTQPVYRVSVPENLPVGTPLAVTATDQDEGVHGEVTS
 FVKITEKISQIFCLNVLTGEISTANLDYEDSSFYELGVEARDGPGLRDRRAKVLITILDVNDNVPEVVVT
 SGSRTIAESAPPGTVIALFQVDRDSSLNGLVTCIPRSLPFELEKSVGNYYRLVTNAALDREEVFLYNI
 TVTATDKGTPPLSTETIISLNVADTNDNPPTFPHSSYSVYVLENNPRGASIFSVNALDPDQNAQVSY
 LAEDTLQGAPLSSYVSINSDTGILYALRSFDYEQRLDLQLWVTASDGDPLSSNVSLSLFVLDQNDNAP
 EILYPALPTDGGTVELAPRSAEPGYLVTKVAVDRDSSGNLAWLSYRLLKASEPGLFSVGLHTGEVTRAR
 ALLDRDALKQSLVVAVQDHGQPPLSATVTLTVAVADRIDILADLGSLEPSAKPNDSLTLYLVVAVAAV
 SCVFLAFVIVLLALRLQRWHKSRLQASGGGLASMPGSHFVGVEGVRAFLQTYSHVSLTADSRKSHLIF
 PQPNYADTLINQESYEKSEPLLITQDLLETKGEPRQLQQAPPNTDWRFSQAQRPQSGSQNGDDTGTWPN
 NQFDTEMLQAMILASASEAADGSSTLGGAGTMGLSARYGPQFTLQHPDYRQNVYIPGSNATLTNAAGK
 RDGKAPAGGNGNKKKSGKKEKK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8010_a03.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:




ACCN: NM_018919


ORF Size: 2796 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_018919.3
RefSeq Size:	4605 bp
RefSeq ORF:	2799 bp
Locus ID:	56109
UniProt ID:	Q9Y5G7
Cytogenetics:	5q31.3
Protein Families:	Transmembrane
MW:	97.7 kDa
Gene Summary:	This gene is a member of the protocadherin gamma gene cluster, one of three related clusters tandemly linked on chromosome five. These gene clusters have an immunoglobulin-like organization, suggesting that a novel mechanism may be involved in their regulation and expression. The gamma gene cluster includes 22 genes divided into 3 subfamilies. Subfamily A contains 12 genes, subfamily B contains 7 genes and 2 pseudogenes, and the more distantly related subfamily C contains 3 genes. The tandem array of 22 large, variable region exons are followed by a constant region, containing 3 exons shared by all genes in the cluster. Each variable region exon encodes the extracellular region, which includes 6 cadherin ectodomains and a transmembrane region. The constant region exons encode the common cytoplasmic region. These neural cadherin-like cell adhesion proteins most likely play a critical role in the establishment and function of specific cell-cell connections in the brain. Alternative splicing has been described for the gamma cluster genes. [provided by RefSeq, Jul 2008]


Product images:



Circular map for RC214360



Western blot validation of overexpression lysate (Cat# [LY412931]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC214360 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified PCDHGA6 protein (Cat# [TP314360]). The protein was produced from HEK293T cells transfected with PCDHGA6 cDNA clone (Cat# RC214360) using MegaTran 2.0 (Cat# [TT210002]).