

Product datasheet for **RC214381L4V**

PHF21A (NM_016621) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PHF21A (NM_016621) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PHF21A
Synonyms:	BHC80; BM-006; IDDBCS; NEDMS
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_016621
ORF Size:	1902 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214381).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_016621.3
RefSeq Size:	7194 bp
RefSeq ORF:	1905 bp
Locus ID:	51317
UniProt ID:	Q96BD5
Cytogenetics:	11p11.2
Domains:	PHD
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)

MW: 70 kDa

Gene Summary: The PHF21A gene encodes BHC80, a component of a BRAF35 (MIM 605535)/histone deacetylase (HDAC; see MIM 601241) complex (BHC) that mediates repression of neuron-specific genes through the cis-regulatory element known as repressor element-1 (RE1) or neural restrictive silencer (NRS) (Hakimi et al., 2002 [PubMed 12032298]).[supplied by OMIM, Nov 2010]