

Product datasheet for **RC214391L2V**

Motilin receptor (MLNR) (NM_001507) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Motilin receptor (MLNR) (NM_001507) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | Motilin receptor |
| Synonyms: | GPR38; MTLR1 |
| Mammalian Cell Selection: | None |
| Vector: | pLenti-C-mGFP (PS100071) |
| Tag: | mGFP |
| ACCN: | NM_001507 |
| ORF Size: | 1236 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC214391). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001507.1 , NP_001498.1 |
| RefSeq Size: | 1239 bp |
| RefSeq ORF: | 1239 bp |
| Locus ID: | 2862 |
| UniProt ID: | O43193 |
| Cytogenetics: | 13q14.2 |
| Protein Families: | Druggable Genome, Transmembrane |
| Protein Pathways: | Neuroactive ligand-receptor interaction |



[View online »](#)

MW: 45.2 kDa

Gene Summary: Motilin is a 22 amino acid peptide hormone expressed throughout the gastrointestinal (GI) tract. The protein encoded by this gene is a motilin receptor which is a member of the G-protein coupled receptor 1 family. This member is a multi-pass transmembrane protein, and is an important therapeutic target for the treatment of hypomotility disorders. [provided by RefSeq, Aug 2011]