

Product datasheet for **RC214431L3V**

CLEC18C (NM_173619) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CLEC18C (NM_173619) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CLEC18C
Synonyms:	MRCL; MRCL3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_173619
ORF Size:	1338 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214431).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_173619.2
RefSeq Size:	2100 bp
RefSeq ORF:	1341 bp
Locus ID:	283971
UniProt ID:	Q8NCF0
Cytogenetics:	16q22.1
MW:	49.6 kDa



[View online »](#)

Gene Summary: Binds polysaccharides in a Ca²⁺-independent manner with a preferential binding to fucoidan, beta-glucans and galactans.[UniProtKB/Swiss-Prot Function]