

Product datasheet for **RC214436A1V**

Human SLAP2 (SLA2) (NM_175077) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human SLAP2 (SLA2) (NM_175077) AAV Particle
Tag: Myc-DDK
Symbol: SLAP2
Synonyms: C20orf156; MARS; SLAP-2; SLAP2
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC214436 representing NM_175077
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGGAAGTCTGCCAGCAGAGAAAATCTCTGCCAAGCCCAAGCTTGAGTTCCTCTGTCCAAGGCCAGG
GACCTGTGACCATGGAAGCAGAGAGAAGCAAGGCCACAGCCGTGGCCCTGGGCAGTTTCCCGGCAGGTGG
CCCGGCCGAGCTGTCGCTGAGACTCGGGGAGCCATTGACCATCGTCTCTGAGGATGGAGACTGGTGGACG
GTGCTGTCTGAAGTCTCAGGCAGAGAGTATAACATCCCCAGCGTCCACGTGGCCAAAGTCTCCCATGGGT
GGCTGTATGAGGGCCTGAGCAGGGAGAAAGCAGAGGAAGTCTGTTGTTACCTGGGAACCTGGAGGGGC
CTTCCTCATCCGGGAGAGCCAGACCAGGAGAGGCTCTTACTCTGTGTCAGTCCGCCTCAGCCGCCCTGCA
TCCTGGGACCGGATCAGACACTACAGGATCCACTGCCTTGACAATGGCTGGCTGTACATCTCACCAGCGCC
TCACCTTCCCCTCACTCCAGGCCCTGGTGGACCATTACTCTGAGGGCTGGCCCGCTCCCTGGCAAGGATA
TACCCCTACCTGTGACTGTGCAGAGGACCACTCAACTGGAAAGAGCTGGACAGCTCCCTCCTGTTTTC

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC214436 representing NM_175077
Red=Cloning site Green=Tags(s)

MGSLPSRRKSLPSPSLSSSVQGGPVTMEAKSKATAVALGSFPAGGPAELSLRLGEPLTIVSEGDWWT
VLSEVSGREYNIPSVHVAKVSHGWL YEGLSREKAEELLLLPGNPGGAFLIRESQTRRGYSLSVRLSRPA
SWDRIRHYRIHCLDNGWLYISPRLTFPSLQALVDHYSEGWPAPWQGYTPTCDCAEDTTQLERAGQLPPVF

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Species: Human



[View online »](#)

Serotype:	AAV-2
ACCN:	NM_175077
ORF Size:	630 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_175077.1</u>
RefSeq Size:	2447 bp
RefSeq ORF:	633 bp
Locus ID:	84174
UniProt ID:	<u>Q9H6Q3</u>
Cytogenetics:	20q11.23
MW:	22.9 kDa