

Product datasheet for **RC214436L3V**

SLAP2 (SLA2) (NM_175077) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | SLAP2 (SLA2) (NM_175077) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | SLAP2 |
| Synonyms: | C20orf156; MARS; SLAP-2; SLAP2 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_175077 |
| ORF Size: | 630 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC214436). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_175077.1 |
| RefSeq Size: | 2447 bp |
| RefSeq ORF: | 633 bp |
| Locus ID: | 84174 |
| UniProt ID: | Q9H6Q3 |
| Cytogenetics: | 20q11.23 |
| Protein Families: | Druggable Genome, Transcription Factors |
| MW: | 22.9 kDa |



[View online »](#)

Gene Summary:

This gene encodes a member of the SLAP family of adapter proteins. The encoded protein may play an important receptor-proximal role in downregulating T and B cell-mediated responses and inhibits antigen receptor-induced calcium mobilization. This protein interacts with Cas-Br-M (murine) ecotropic retroviral transforming sequence c. Two transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Jul 2008]