

## Product datasheet for RC214445

### STAG3L4 (NM\_022906) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** STAG3L4 (NM\_022906) Human Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** STAG3L4  
**Synonyms:** FLJ13195  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >RC214445 representing NM\_022906  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGAGTGGCTGGATAGCTACAAGCAAGACCAGGATGCAGGATTTCTGGAGCTTGTTAACTTTTTCATCTG  
ATCTTGTGGATGTAAGGACTCGGGGACTATCCTCTGACAGCTCCAGGTCTATCCTGGAAGAAGTTCCA  
GGGCAGCTTCTGTGAGTTTGTGGGACATTGGTCTGTCGGTGCCAGTACATCCTCCTCCATGATGACTTC  
CCTATGGACAACCTCATCTCCCTGCTCACTGGCTTCTCAGACTCACAAGTCTGCGCTTCTGTCACTA  
GCACCCTGGCTGCTATGAACTGATGACCTCCCTGGTAAGAGTTGCCCTCCAAGTCTGCACCAAGA  
TATCAATCAGCGTCAGTATGAGGCTGAAAGAAACAAGGGGCCAGGGCAGAGGGCACCTGAGCGGCTGGAG  
AGCCTGTTGGAGAAACACAAAGAGCTGCAT

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC214445 representing NM\_022906  
Red=Cloning site Green=Tags(s)

MSGWIATSKTRMQDFWSLLTFSSDLVDVKDSGDYPLTAPGLSWKKFQGSFCEVGTLVCRQYILLHDDF  
PMDNLISLLTGFSDSQVCAFCSTLAAMKLMTSLVRVALQLSLHQDINQRQYEAERNKGPGRAPERLE  
SLLEKHKELH

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/mk6486\\_e09.zip](https://cdn.origene.com/chromatograms/mk6486_e09.zip)



[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM\_022906

ORF Size: 452 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_022906.1](#), [NP\\_075057.1](#)


RefSeq Size: 2147 bp

RefSeq ORF: 452 bp


Locus ID: 64940

Cytogenetics: 7q11.21  
 MW: 16.8 kDa

**Product images:**



Circular map for RC214445



Western blot validation of overexpression lysate (Cat# [LY411479]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC214445 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).