

## Product datasheet for **RC214463L3V**

### LSS (NM\_001001438) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	LSS (NM_001001438) Human Tagged ORF Clone Lentiviral Particle
Symbol:	LSS
Synonyms:	APMR4; CTRCT44; HYPT14; OSC
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001001438
ORF Size:	2196 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC214463).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001001438.1</a> , <a href="#">NP_001001438.1</a>
RefSeq Size:	2658 bp
RefSeq ORF:	2199 bp
Locus ID:	4047
UniProt ID:	<a href="#">P48449</a>
Cytogenetics:	21q22.3
Protein Families:	Druggable Genome
Protein Pathways:	Metabolic pathways, Steroid biosynthesis



[View online »](#)

**MW:** 83.1 kDa

**Gene Summary:** The protein encoded by this gene catalyzes the conversion of (S)-2,3 oxidosqualene to lanosterol. The encoded protein is a member of the terpene cyclase/mutase family and catalyzes the first step in the biosynthesis of cholesterol, steroid hormones, and vitamin D. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Feb 2009]