

Product datasheet for **RC214464**

IL4R (NM_000418) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	IL4R
Synonyms:	CD124; IL-4RA; IL4RA
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>RC214464 representing NM_000418.

Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAAATTTTGAATACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGGGTGGCTTTGCTCTGGGCTCCTGTTCCCTGTGAGCTGCCTGGTCTGCTGCAGGTGGCAAGCTCT
GGGAACATGAAGTCTTGCAGGAGCCACCTGCGTCTCCGACTACATGAGCATCTCTACTTGCAGTGG
AAGATGAATGGTCCACCAATTGCAGCACCGAGCTCCGCCTGTTGTACCAGCTGGTTTTTCTGCTCTCC
GAAGCCACACGTGTATCCCTGAGAACAACGGAGGCGCGGGGTGCGTGTGCCACCTGCTCATGGATGAC
GTGGTCACTGCGGATAACTATACACTGGACCTGTGGGCTGGGCAGCAGCTGCTGTGAAGGGCTCCTTC
AAGCCAGCGAGCATGTGAAACCCAGGGCCCCAGGAACTGACAGTTACACCAATGTCTCCGACACT
CTGCTGCTGACCTGGAGCAACCCGTATCCCCCTGACAATTACCTGTATAATCATCTACCTATGCAGTC
AACATTTGGAGTGAACGACCCGGCAGATTTTCAAGTCTATAACGTGACCTACCTAGAACCTCCCTC
CGCATCGCAGCCAGCACCTGAAGTCTGGGATTTCTACAGGGCACGGGTGAGGGCTGGGCTCAGTGC
TATAACACCACCTGGAGTGAAGTGGAGCCCCAGCACCAAGTGGCACAACCTCTACAGGGAGCCCTTCGAG
CAGCACCTCCTGCTGGGCTCAGCGTTTCTGCAATTGTATCCTGGCCGCTGCTGTTGTGCTATGTC
AGCATACCAAGATTAAGAAAGAAATGGTGGGATCAGATTTCCCAACCCAGCCCGCAGCCGCTCGTGGCT
ATAATAATCCAGGATGCTCAGGGGTACAGTGGGAGAAGCGGTCCCGAGGCCAGGAACAGCCAAGTGC
CCACACTGGAAGAATTGTCTTACCAAGCTCTTGCCCTGTTTTCTGGAGCACAACATGAAAAGGGATGAA
GATCCTCACAAGGCTGCCAAAGAGATGCCTTTCCAGGGCTCTGAAAAATCAGCATGGTGGCCAGTGGAG
ATCAGCAAGACAGTCTCTGGCCAGAGAGCATCAGCGTGGTGGATGTGTGGAGTTGTTTGGGCCCCG
GTGGAGTGTGAGGAGGAGGAGGTAGAGGAAGAAAAAGGGAGCTTCTGTGCATCGCCTGAGAGCAGC
AGGGATGACTTCCAGGAGGAAGGGAGGGCATTGTGGCCCGGCTAACAGAGAGCCTGTTCTGGACCTG
CTCGGAGAGGAGAATGGGGGCTTTTGCCAGCAGGACATGGGGGAGTCATGCCTTCTCCACCTTCGGGA
AGTACGAGTGCTCACATGCCCTGGGATGAGTTCCTCAAGTGCAGGGCCCAAGGAGGCACCTCCCTGGGGC
AAGGAGCAGCCTCTCCACCTGGAGCCAAGTCTCTGCCAGCCCGACCCAGAGTCCAGACAACCTGACT
TGCACAGAGACGCCCCCTCGTCATCGCAGGCAACCTGCTTACCGCAGTTTACGCAACTCCCTGAGCCAG
TCACCGTGTCCCAGAGAGCTGGGTCCAGACCCACTGCTGGCCAGACACCTGGAGGAAGTAGAACCCGAG
ATGCCCTGTGTCCCCAGCTCTCTGAGCCAACCACTGTGCCCCAACCTGAGCCAGAAACCTGGGAGCAG
ATCCTCCGCCGAAATGTCTCCAGCATGGGGCAGCTGCAGCCCCGTCTCGGCCCCACCACTGGGTAT
CAGGAGTTTGTACATGCGGTGGAGCAGGTGGCACCCAGGCCAGTGGGTGGTGGCTTGGGTCCCCCA
GGAGAGGCTGGTTACAAGGCCTTCTCAAGCCTGCTTGCCAGCAGTGTGTGTCCCAGAGAAATGTGGG
TTTGGGGCTAGCAGTGGGGAAGAGGGGTATAAGCCTTTCCAAGACCTCATTCTGGCTGCCTGGGGAC
CCTGCCCCAGTCCCTGTCCCCTTGTTCACCTTTGGACTGGACAGGGAGCCACCTCGCAGTCCGACAGC
TCACATCTCCCAAGCAGTCCCCAGAGCACCTGGGTCTGGAGCCGGGGGAAAAGGTAGAGGACATGCCA
AAGCCCCCACTTCCCCAGGAGCAGGCCACAGACCCCTTGTGGACAGCCTGGGCAGTGGCATTGTCTAC
TCAGCCCTTACCTGCCACCTGTGCGGCCACCTGAAACAGTGTATGGCCAGGAGGATGGTGGCCAGACC
CCTGTCTATGGCCAGTCTTGTGTGGTGTCTGTGTGGAGACAGGTCTCGCCCCCTACAACCCCTG
AGGGCCCCAGACCCCTCTCCAGGTGGGGTTCCACTGGAGGCCAGTCTGTGTCCGGCCTCCCTGGCACCC
TCGGGCATCTCAGAGAAGAGTAAATCCTCATCATCCTTCCATCCTGCCCTGGCAATGCTCAGAGCTCA
AGCCAGACCCCCAAAATCGTGAACCTTTGTCTCCGTGGGACCCACATACATGAGGGTCTCT
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Protein Sequence: >Peptide sequence encoded by RC214464
Blue=ORF Red=Cloning site Green=Tag(s)



MGWLCSGLLPVSVCLVLLQVASSGNMKVLQEPTCVSDYMSISTCEWMNGPTNCSTELRLLYQLVFLLS
EAHTCIPENGGAGCVCHLLMDDVVSADNYTLDLWAGQQLLWGSFKPSEHVKPRAPGNLTVHTNVSDT
LLL TWSNPYPDPNYL YNHL TYAVNIWSENDPADFRIYNVTYLEPSLRIAASTLKSIGSYRARVRAWAQC
YNTTWSEWSPSTKWHNSYREPFEQHLLLVGSVSCIVILAVCLLYCVYSITKIKKEWWDQIPNPARSRLVA
II IQDAQGSQWEKRSRGQEPAKCPHWKNCLETKLLPCFLEHNMKRDEDPHAKAEMPFGSGKSAWCVPVE
ISKTVLWPESISVVRVCVELFEAPVECEEEEEVEEEKGSFCASPESSRDDFQEGREGIVARL TESLFLDL
LGEENGFCQQDMGESCLLPSPSGTSAHMPWDEFPSAGPKEAPPWGKEQLHLEPSPASPTQSPDNL
CTETPLVIAGNPAYRSFSNLSQSPCPRELGPDPLLARHLEEEPEMPCVPQLSEPTTVPQPEPETWEQ
ILRRNVLQHGA AAPVSAPTSYG YEFVHAEQGGTQASAVVGLGPPGEAGYAFSSLLASSAVSPEKCG
FGASSGEEGYKPFQDLIPGCGDPAPVPVPLFTFGLDREPPRSPQSSHLPSSSPEHLGLEPGEKVEDMP
KPPLPQEQTADPLVDSLGSIGIYSSALTEKLCGHLLKQCHGQEDGGQTPVMASPCCGCCCGDRSSPPTPL
RADPSPGGVPLEASLCPASLAPSGTSEKSKSSSFHPAPGNAQSSSQTPKIVNFVSVGPTYMRVS
TRTRPLEOKLISEEDLAANDILDYKDDDDKY

Chromatograms: https://cdn.origene.com/chromatograms/mk6237_f04.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM 000418


ORF Size: 2475 bp

OTI Disclaimer:	Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery. The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_000418.4
RefSeq Size:	3678 bp
RefSeq ORF:	2478 bp
Locus ID:	3566
UniProt ID:	P24394
Cytogenetics:	16p12.1
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	Cytokine-cytokine receptor interaction, Hematopoietic cell lineage, Jak-STAT signaling pathway
MW:	89.7 kDa


Gene Summary:

This gene encodes the alpha chain of the interleukin-4 receptor, a type I transmembrane protein that can bind interleukin 4 and interleukin 13 to regulate IgE production. The encoded protein also can bind interleukin 4 to promote differentiation of Th2 cells. A soluble form of the encoded protein can be produced by proteolysis of the membrane-bound protein, and this soluble form can inhibit IL4-mediated cell proliferation and IL5 upregulation by T-cells. Allelic variations in this gene have been associated with atopy, a condition that can manifest itself as allergic rhinitis, sinusitis, asthma, or eczema. Polymorphisms in this gene are also associated with resistance to human immunodeficiency virus type-1 infection. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Apr 2012]

Product images:



Circular map for RC214464



Western blot validation of overexpression lysate (Cat# [LY424733]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC214464 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).