

Product datasheet for **RC214467A1V**

Human FCRL3 (NM_052939) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human FCRL3 (NM_052939) AAV Particle
Tag:	Myc-DDK
Symbol:	FCRL3
Synonyms:	CD307c; FCRH3; IFGP3; IRTA3; SPAP2
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)



[View online »](#)

ORF Nucleotide Sequence:

```
>RC214467 representing NM_052939
Red=Cloning site Blue=ORF Green=Tags(s)
```

CTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

[illegible]

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCTGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:	>RC214467 representing NM_052939 Red=Cloning site Green=Tags(s) MLLWLLLLILTPGREQSGVAPKAVLLLNPWPSTAFKGEKVALICSSISHSLAQGDYWYHDEKLLKIKHD KIQITEPGNYQCKTRGSSLSDAVHVEFSPDWLILQALHPVFEGDNVILRCQKDNKNTHQKVVYKDGKQL PNSYNLEKITVNSVSRDNSKYHCTAYRKFYILDIEVTSKPLNIQVQELFLHPVLRASSSTPIEGSPMTLT CETQLSPQRPDVQLQFSLFRDSQTLGLGWSRSPRLQIPAMWTEGSGSYWCEVETVTHSIKKRSLRSQIRV QRVPVSNVNLEIRPTGGQLIEGENMVLICSVAQSGTVTFSWHKEGRVRSGLGRKTQRSLLAELHVLTVKE SDAGRYYCAADNVHSPILSTWIRVTVRIPVSHPLVTFRAPRAHTVVGDLLELHCESLRGSPPILYRFYHE DVTLGNSSAPSGGGASFNLSTAETHSGNYSCDADNGLGAQHSVGSLRVTPVSRPVLTLRAPGAQAVVG DLLELHCESLRGSPPILYWYFHEDDTLGNISAHSGGGASFNLSTTEHSGNYSCEADNGLGAQHSKVVTL NVTGTSRNRGTGLTAAGITGLVLSILVLAALLHYARARRKPGGLSATGTSSHSPSECQEPSSSRPSRI DPQEPHSLKPLAPMELEPMYSNVNPGDSNPIYSQIWSIQHTKENSANCPMMHQEHEELTVLYSELKKTHTP DDSAGEASSRGRAHEEDDEENYENVPRVLLASDH TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_052939
ORF Size:	2202 bp
Buffer:	PBS with 0.001% Pluronic F7219
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_052939.3
RefSeq Size:	3019 bp
RefSeq ORF:	2205 bp
Locus ID:	115352
UniProt ID:	Q96P31
Cytogenetics:	1q23.1
MW:	79.6 kDa